



NJW SERIES

High Precision Balance
0.01g~0.1g Display Resolution





FEATURES

- 1/30000 Display Resolution High Precision Balance
- Suitable for Laboratory and Jewelry Application
- 13 User Programmable Weight Units: - Metric, Metric Carat, Pound (Avoirdupois), Ounce (Avoirdupois), Dram (Avoirdupois), Grain (U.K.), Troy Ounce, Pennyweight, Momme (Japan), Tael (Taiwan), Tael (China), Troy Tael, Tola (India)
- Equipped with Power On Zero Tracking, Manual Zero Setting, Power Saving, Auto Calibration, Percentage and Counting Function
- Long Battery Operating Time: 100 Hours Plus
- Dual Power: By 4 X UM-2 Batteries or Power Adaptor
- Low Battery Warning Signal
- 6 x 13 mm Bold Type Wide Angle Display Digits
- Power Adaptor, Backlight, Adjustable Feet and Dust Cover Included
- NJW-300 and NJW-600 Equipped with Bubble Level and Wind Shield

OPTION AVAILABLE

- RS232C Interface Output

SPECIFICATIONS

Model No.	Capacity	Division	Platter Size
NJW-300	300g	0.01g	ø110mm
NJW-600	600g	0.02g	ø110mm
NJW-1500	1500g	0.05g	140 x 170mm
NJW-3000	3000g	0.1g	140 x 170mm
Power Source	By 4 x UM-2 Batteries or Power Adaptor		
Operation Environment	0~40°C (32~104°F), Non-condensed R.H. ≤ 85%		
Individual Package	NJW-300 & 600: 310(L) x 300(W) x 170(H)mm Other Models: 360(L) x 275 x(W) 100(H)mm		
Net/Gross Weight	NJW-300 & 600: 1.22/2.00kg Other Models: 1.24/1.60kg		

• Specifications subject to change without notice



UNIVERSAL WEIGHT ENTERPRISE CO., LTD.

2-5 Floor No. 39 Pao Shing Road Hsin Tien City 231 Taipei Hsien Taiwan R.O.C.
Tel : +886-2-29180121 (12 Lines) Fax : +886-2-29181508 · 29183652 · 29112176-7
http://www.uwe-scale.com.tw e-mail:uwe@msb.hinet.net

Apr '99 Designed by DAILY
Printed in Taiwan